

Total Alkalinity

SEAL Application ML-905

SCOPE AND APPLICATION

This method covers the determination of alkalinity in wastewater, surface and drinking waters.

Applicable Range: 20 -


20 - 1000 mg/L as $CaCO_3$

This method is equivalent to Standard Method 2320 B, ASTM D 1067 and ISO 9963-2.

METHOD DESCRIPTION

Hydroxyl ions present in a sample as a result of dissociation or hydrolysis of solutes react with additions of standard acid. Alkalinity determined based on the end-point pH used.

METHOD PERFORMANCE



The titration solution used was 0.1 N hydrochloric acid.

Alkalinity Method Performance mg/L CaCO₃

g, 2 00003		
20.7	201.7	493.7
19.3	198.9	495.6
20.0	198.6	496.5
20.5	201.5	499.1
20.5	200.6	493.7
19.3	197.2	504.6
19.2	196.3	511.7
19.4	201.0	496.9
19.3	193.0	497.3
True Value		
20.0	200.0	500.0
Average		
19.7	198.8	498.8
*Standard Deviation		
0.5	2.7	5.5
**Relative Standard Deviation		
2.8%	1.4%	1.1%
Recovery		
98.8%	99.4%	99.8%

^{*} Standard deviation based on the entire population.

^{**} Relative standard deviation calculated by dividing the standard deviation by the mean and multiplying by 100.